




Curve "0" SO1 Before 12/1970 SEV-Marchal 41 101 012 C114/1 2,7L  
 Check Timing: 2000 RPM / 27 Degrees Static Timing: 800 RPM / 1 Degree



Curve "1" SM SB 12/1970 → 3/1971 SEV-Marchal 41 101 012 C114/1 2,7L  
 Check Timing: 2000 RPM / 27 Degrees Static Timing: 900 RPM / 4 Degrees




Curve "2" SO2 Merak/Merak SS SEV-MARCHAL 411 10 902 C114/12 3.0L Euro  
 Check Timing: 2200 RPM / 22 Degrees Static Timing: 700 RPM / 2 Degrees



Curve "3"      SO2 & RD2 Vacuum Retard      SEV-Marchal 411 10 902      C114/11 3,0L Automatic  
 C114/13 3,0L De-Poluted

Check Timing: 2000 RPM / 26 Degrees      Static Timing: 700 RPM / 6 Degrees



Curve "4"      SU1      C114/1 2,7L De-Poluted

Check Timing: 2000 RPM / 27 Degrees      Static Timing: 900 RPM / 1 Degree



Curve "5"      SI1 & RD1 Vacuum Retard      SEV-Marchal 4114 0004      C114/03 2,7L Injection

Check Timing: 2000 RPM / 22 Degrees      Static Timing: 925 RPM / 10 Degrees

Curve selection:



Remark:

Although this information was composed with the greatest possible diligence, Albertronic BV cannot accept any liability for consequences arising from the use of this information.